

# BED SETTEE MULTY

## Designer

Claude Brisson



## Dimensions

## Technical characteristics

Frame in sambawood, multi-ply and panels of particleboard, with steel mechanism and birch slats. Seat in polyether foam ( $28\text{kg/m}^3$  -  $3.0\text{kPa}$ ). Back in polyether foam ( $28\text{kg/m}^3$  -  $3.0\text{kPa}$ ). Lumbar cushions are supplied as follows: 2 x W60 for small bed-settee without arms, and 2 x W78 for large bed-settee without arms; 1 x W78 for small bed-settee with arms, and 2 x W60 for large bed-settee with arms. Please take these details into consideration when ordering additional/subsequent lumbar. The armrests (polyether foam  $28\text{ kg/m}^3$  -  $3.0\text{ kPa}$  + polyester quilting  $200\text{ g/m}^2$ ) or bedside tables (MDF top diameter  $43\text{ cm}$  - thickness  $2\text{ cm}$ ) pivot to allow the bed-settee to open.

MULTY is easy to operate, and offers excellent comfort both as a bed (with polyether mattress size :  $153\times 194\times 8$ ,  $119\times 194\times 8$ ) and as a settee. When unfolded, its depth is  $200\text{ cm}$  for the bed-settee W  $155$ . Mattress has an ecru-coloured cover in  $50\%$  polyester/ $50\%$  cotton (dry cleanable), which has been treated with "HealthGuard" (antibacterial & fungicidal). A foam sheet is supplied to place between the slats and the seat when folded. For instructions on use, please see the card which comes with the bed. MULTY has been specially designed for use as a bed on an everyday basis. Optional extras: lumbar cushions and arms. Covers of lumbar cushions and armrests are removable.